



Recombinant Human TNFAIP3-interacting protein 1 (TNIP1)

Product Code	CSB-MP619959HU
Abbreviation	TNIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15025
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEGRGPYRIY DPGGSVPSGE ASAAFERLVK ENSRLKEKMQ GIKMLGELLE ESQMEATRLR QKAEELVKDN ELLPPPSPSL GSFDP LAELT GKDSNVTASP TAPACPSDKP APVQKPPSSG TSSEFEVVTP EEQNSPESSS HANAMALGPL PREDGNLMLH LQRLETTLSV CAEEDPHGQL FTHLGRMALE FNRLASKVHK NEQRTSILQT LCEQLRKENE ALKAKLDKGL EQRDQAAERL REENLELKKL LMSNGNKEGA SGRPGSPKME GTGKKAVAGQ QQASVTAGKV PEVVALGAAE KVKMLEQQR SELLEVNKQW DQHFRSMKQQ YEQKITELRQ KLADLQKQVT DLEAEREQKQ RDFDRKLLLA KSKIEMEETD KEQLTAEAKE LRQKVYQLQD QLSPLTRQRE YQEKEIQRLN KALEEALSIQ TPPSSPPTAF GSPEGAGALL RKQELVTQNE LLKQQVKIFE EDFQRERSDR ERMNEEKEEL KKQVEKLQAQ VTLSNAQLKA FKDEEKAREA LRQQKRKAKA SGERYHVEPH PEHLGAYPY AYPPMPAMVP HHGFEDWSQI RYPPPPMAME HPPPLPNSRL FHLPEYTWRL PCGGVRNPQ SSQVMDPPTA RPTEPESPKN DREGPQ
Source	Mammalian cell
Target Names	TNIP1
Protein Names	Recommended name: TNFAIP3-interacting protein 1 Alternative name(s): HIV-1 Nef-interacting protein Nef-associated factor 1 Short name= Naf1 Nip40-1 Virion-associated nuclear shuttling protein Short name= VAN Short nam
Expression Region	1-636
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.