



Recombinant Human Histone-lysine N-methyltransferase SETDB1 (SETDB1)

| | |
|--------------------------|--|
| Product Code | CSB-EP619960HU-B |
| Abbreviation | SETDB1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q15047 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MSSLPGCIGLDAATATVESEEIAELQQAVVEELGISMEELRHFIDEELEKMDCV QQRKKQLAELETWVIQKESEVAHVVDQLFDDASRAVTNCESLVKDFYSKLGLQ YRDSSSEDESSRPTEIIEIPDEDDDDVLSIDSGDAGSRTPKDQKLREAMAALRKS AQDVQKFMDAVNKKSSSQDLHKGTLSQMSGELSKDGDIVSMRILGKKRKT WHKGTLIAIQTVGPGKKYKVKFDNKGKSLLSGNHAIYDYHPPADKLYVGSRVV AKYKDGNQVWLYAGIVAETPNVKNKLRFLIFFDDGYASYVTQSELYPICRPLKK TWEDIEDISCRDFIEEYVTAYPNRPMVLLKSGQLIKTEWEGTWWKSRVEEVDG SLVRILFLVLFSTILEAEVGGGGT |
| Source | E.coli |
| Target Names | SETDB1 |
| Protein Names | Recommended name: Histone-lysine N-methyltransferase SETDB1 EC= 2.1.1.43 Alternative name(s): ERG-associated protein with SET domain Short name= ESET Histone H3-K9 methyltransferase 4 Short name= H3-K9-HMTase 4 Lysine |
| Expression Region | 1-397 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length of Isoform 2 |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.