



# Recombinant Human Vacuolar protein sorting-associated protein 72 homolog (VPS72)

<b>Product Code</b>	CSB-MP620899HU
<b>Abbreviation</b>	VPS72
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15906
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSLAGGRAPR KTAGNRLSGL LEAEEDEFY QTTYGGFTEE SGDDEYQGDO SDTEDEVSD FDIDEGDEPS SDGEAEEPRR KRRVVTKAYK EPLKSLRPRK VNTPAGSSQK AREEKALLPL ELQDDGSDSR KSMRQSTAEH TRQTFRLVQE RQGQSRRRKG PHCERPLTQE ELLREAKITE ELNLRSLITY ERLEADKKKQ VHKRKCPCGP IITYHSVTVP LVGEPGPKKE NVDIEGLDPA PSVSALTPHA GTGPNVPPAR CSRTFITFSD DATFEEWFPQ GRPPKVPVRE VCPVTHRPAL YRDPVTDIPY ATARAFKIIR EAYKKYITAH GLPPTASALG PGPPPPEPLP GSGPRALRQK IVIK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	VPS72
<b>Protein Names</b>	Recommended name: Vacuolar protein sorting-associated protein 72 homolog Alternative name(s): Protein YL-1 Transcription factor-like 1
<b>Expression Region</b>	1-364
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.