



Recombinant Human Insulin-like growth factor-binding protein 7 (IGFBP7)

Product Code	CSB-BP620956HU
Abbreviation	IGFBP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16270
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SSSD TCGPCEPASC PPLPPLGCLL GETRDACGCC PMCARGECEP CGGGGAGRGY CAPGMECVKS RKRRKKGAGA AAGGPGVSGV CVCKSRYPVC GSDGTTYPSG CQLRAASQRA ESRGEKAITQ VSKGTCEQGP SIVTPPKDIW NVTGAQVYLS CEVIGIPTPV LIWNKVKRGH YGVQRTELLP GDRDNLAIQT RGGPEKHEVT GWVLVSPLSK EDAGEYECHA SNSQQQASAS AKITVVDALH EIPVKKGEGA EL
Source	Baculovirus
Target Names	IGFBP7
Protein Names	Recommended name: Insulin-like growth factor-binding protein 7 Short name= IBP-7 Short name= IGF-binding protein 7 Short name= IGFBP-7 Alternative name(s): IGFBP-rP1 MAC25 protein PGI2-stimulating factor Prostacycl
Expression Region	27-282
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.