



# Recombinant Bovine AP-2 complex subunit sigma (AP2S1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP621049BO                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Abbreviation</b>      | AP2S1                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                             |
| <b>Uniprot No.</b>       | Q17QC5                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | Bos taurus (Bovine)                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MIRFILIQNR AGKTRLAKWY MQFDDDEKQK LIEEVHAVVT VRDAKHTNFV<br>EFRNFKIIYR RYAGLYFCIC VDVNDNNLAY LEAIHNFVEV LNEYFHNVCE<br>LDLVFNFYKV YTVVDEMFLA GEIRETSQTK VLKQLLMLQS LE                                                                                                                                                                                                                                          |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Target Names</b>      | AP2S1                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: AP-2 complex subunit sigma Alternative name(s): Adapter-related protein complex 2 sigma subunit Adaptor protein complex AP-2 subunit sigma Clathrin assembly protein 2 small chain Clathrin coat assembly protein AP17                                                                                                                                                                    |
| <b>Expression Region</b> | 1-142                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Target Details</b>    | One of two major clathrin-associated adaptor complexes, AP-2, is a heterotetramer which is associated with the plasma membrane. This complex is composed of two large chains, a medium chain, and a small chain. This gene encodes the small chain of this complex. Alternative splicing has been observed in this gene and results in two known transcripts.                                               |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.                                                                                                                                                                                                                                                                  |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.