



Recombinant F-box/SPRY domain-containing protein 1 (fsn-1)

Product Code	CSB-MP621106CXY
Abbreviation	fsn-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q18223
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAENDGETIV PDEQCNLTA S TPMKSSGLEQ VESPKERKTS MVECDDEGNP SSGTPDESPK QLTPHSIPPR RRRSPRRPEV SASRLPLKVL NQIFQYLPLK DLRSAMLTCH SWNNALSMED SDIWQYLLGK KLPEAAVSDP FLLAELGSAK KKLRAWYFAW NTSDISRNNY IRTNGFTVHR QPVAQSTDGV RGKRGISKGV HAFDITWDGP LGTVAVVGIA TKHAALHCVG YVALLGSDDQ SWGWNLVDNV LMHNGAQLGV YPKMNNPPKY EVGDKIRLII DCDTHVAYFE RNSEFLGIAF NHIPPLRLYP AVCAVYGNT E VTMVYVGSPQ MG
Source	Mammalian cell
Target Names	fsn-1
Protein Names	Recommended name: F-box/SPRY domain-containing protein 1 Alternative name(s): F-box synaptic protein 1
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.