



Recombinant Putative pyridoxal-phosphate dependent protein F13B12.4 (F13B12.4)

Product Code	CSB-BP621177CXY
Abbreviation	F13B12.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19374
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	ED ITKISESRIP NEDAVPDYEA KWRLGAIQKL WNERRKMGHT PMTKFSPPGF PNADIFFKNE TATATRTLKH RFAWALLLWA ITEGKVTSKT SAVYDSTSGN TGSAEAYMCT LVNVPYYAVV ADNLEKEKVK QIESFGGKII KVPVALRNAK AKEFADKNHG FYMNQFGNAE KAEEFHESGD FYFESTNVYH EIIVQLKKDK EQIVKIPDYF VHSAGTGGTI SSVGRYVARY GAPTKVVLSLSD SQYSLFYDYV IGHKFTNQSG AGIWTPPGIA GIGYGYDIEP VWYGETTSLT RNVIHEAMKM PDIASVATMR ILDEKGYNVG PSTSLNFLVS LYKAYQNKAR KSAIKHRLTI VTLACDPGDF YRSTYLNNEW VEKSFKKFGG VMGMECWKKL IQESIDTGSD FYSKGLTMCP GAFKV
Source	Baculovirus
Target Names	F13B12.4
Protein Names	Recommended name: Putative pyridoxal-phosphate dependent protein F13B12.4
Expression Region	19-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.