



Recombinant *Saccharomyces cerevisiae*

Ribosome biogenesis protein SSF2 (SSF2)

Product Code	CSB-BP621459SVG
Abbreviation	SSF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12153
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRRQKKRT HAQITPEQER DIPKSMVIRV GQTSLANHSL NQLVKDFRQI MQPHTAVKLLK ERKSNKLDKF VVMCGPLGVT HLFMFTQSEK TGNVSLKIAR TPQGPTVTFQ VLDYSLGRDI KKFLKRPKSL NNDDVLNPPL LVLNGFSTSK RSEDDQDVN VEKVIVSMFQ NIFPPLNPAR TSLNSIKRIF MINKDRETGE ISMRHYFIDI REVEISRNLK RLYKAKNNLS KTVPNLHRKE DISSLILDHD LGAYTSESEI EDDAIVRVVD NQDVKAKHSQ TSLSQKTPVK MTDNEEREKG IEEEDVEMEE PKPSENSQPT PRKKAIKLTE LGPRLTLKLV KIEDGICSGK VLHHEFVQKS SEEIKALEKR HAAKMRLKEQ RRKEQEENIA KKKAVKDAKK QRKLERRKAR AEEQGEGQGK DGAMSDDGSS SSEDEHYSDV PEDLDSDLFS EVE
Source	Baculovirus
Target Names	SSF2
Protein Names	Recommended name: Ribosome biogenesis protein SSF2
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.