



Recombinant Human MDS1 and EVI1 complex locus protein MDS1 (MECOM)

Product Code	CSB-YP621671HU
Abbreviation	MECOM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13465
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRSKGRARKL ATNNECVYGN YPEIPLEEMP DADGVASTPS LNIQEPCSPA TSSEAFTPKE GSPYKAPIYI PDDIPIPAEF ELRESNMPGA GLGIWTKRKI EVGEKFGPYV GEQRSNLKDP SYGWEVHLPR SRRVSVHSWL YLGKRSSDVG IAFSQADVYM PGLQCAFLS
Source	Yeast
Target Names	MECOM
Protein Names	Recommended name: MDS1 and EVI1 complex locus protein MDS1 Alternative name(s): Myelodysplasia syndrome 1 protein Myelodysplasia syndrome-associated protein 1
Expression Region	1-169
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.