



# Recombinant Human snRNA-activating protein complex subunit 2 (SNAPC2)

<b>Product Code</b>	CSB-EP621673HU-B
<b>Abbreviation</b>	SNAPC2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13487
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKPPPRRRAA PARYLGEVTG PATWSAREKR QLVRLQARQ GQPEPDATEL ARELRGRSEA EIRVFLQQLK GRVAREAIQK VHPGGLQGPR RREAQPPAPI EVWTDLAEKI TGPLEEALAV AFSQVLTIAA TEPVTLHLSK PPKPTQARGK PLLSAPGGQ EDPAPEIPSS APAAPSSAPR TDPAPEKPS ESSAGPSTEE DFAVDFEKIY KYLSSVSRSG RPELSAAES AVVLDLLMSL PEELPLLCT ALVEHMTETY LRLTAPQPIP AGGSLGPAAE GDGAGSKAPE ETTPATEKAE HSELKSPWQA AGICPLNPFL VPLELLGRAA TPAR
<b>Source</b>	E.coli
<b>Target Names</b>	SNAPC2
<b>Protein Names</b>	Recommended name: snRNA-activating protein complex subunit 2 Short name= SNAPc subunit 2 Alternative name(s): Proximal sequence element-binding transcription factor subunit delta Short name= PSE-binding factor subunit delta Short
<b>Expression Region</b>	1-334
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a subunit of the snRNA-activating protein complex which is associated with the TATA box-binding protein. The encoded protein is necessary for RNA polymerase II and III dependent small-nuclear RNA gene transcription. Alternative splicing results in multiple transcript variants.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.