



Recombinant Human Caspase-8 (CASP8), partial

Product Code	CSB-EP621791HU-B
Abbreviation	CASP8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14790
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SESQ TLDKVVYQMK S KPRGYCLIIN NHNFAKAREK VPKLHSIRDR NGTHLDAGAL TTTFEELHFE IKPHDDCTVE QIYEILKIYQ LMDHSNMDCF ICCILSHGDK GIIYGTGQEQ APIYELTSQF TGLKCPFLAG KPKVFFIQAC QGDNYQKGIP VETD
Source	E.coli
Target Names	CASP8
Protein Names	Recommended name: Caspase-8 Short name= CASP-8 EC= 3.4.22.61 Alternative name(s): Apoptotic cysteine protease Apoptotic protease Mch-5 CAP4 FADD-homologous ICE/ced-3-like protease FADD-like ICE Short name= FLICE ICE-I
Expression Region	217-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.