



Recombinant Human Nicotinate-nucleotide pyrophosphorylase [carboxylating] (QPRT)

Product Code	CSB-EP621868HU-B
Abbreviation	QPRT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15274
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDA EGLALLL PPVTLAALVD SWLREDCPGL NYAALVSGAG PSQAALWAKS PGVLAGQPFF DAIFTQLNCQ VSWFLPEGSK LVPVARVAEV RGP AHCLLLG ERVALNTLAR CSGIASAAAA AVEAARGAGW TGHVAGTRKT TPGFRLVEKY GLLVGGAASH RYDLGGLVMV KDNHVVAAGG VEKAVRAARQ AADFTLKVEV ECSSLQEAVQ AAEAGADLVL LDNFKPEELH PTATVLKAQF PSVAVEASGG ITLDNLPQFC GPHIDVISMG MLTQAAPALD FSLKLFKEV APVPKIH
Source	E.coli
Target Names	QPRT
Protein Names	Recommended name: Nicotinate-nucleotide pyrophosphorylase [carboxylating] EC= 2.4.2.19 Alternative name(s): Quinolate phosphoribosyltransferase [decarboxylating] Short name= QAPRTase Short name= QPRTase
Expression Region	1-297
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a key enzyme in catabolism of quinolate, an intermediate in the tryptophan-nicotinamide adenine dinucleotide pathway. Quinolate acts as a most potent endogenous exitotoxin to neurons. Elevation of quinolate levels in the brain has been linked to the pathogenesis of neurodegenerative disorders such as epilepsy, Alzheimer s disease, and Huntington s disease.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.