



Recombinant Human Protein SSXT (SS18)

Product Code	CSB-EP621876HU
Abbreviation	SS18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15532
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SVAFAAPRQ RGKGEITPAA IQKMLDDNNH LIQCIMDSQN KGKTSECSQY QQMLHTNLVY LATIADSNQN MQSLLPAPPT QNMPMGPGGM NQSGPPPPPR SHNMPSDGMV GGGPPAPHMV NQMNGQMPGP NHMPMQGPGP NQLNMTNSSM NMPSSSHGSM GGYNHSVPSS QSMPVQNQMT MSQGQPMGNY GPRPNMSMQP NQGPMMHQQP PSQQYNMPQG GGQHYQGQQP PMGMMGQVNV GNHMMGQRQI PPYRPPQQGP PQQYSGQEDY YGDQYSHGGQ GPPEGMNQQY YPDGHNDYGY QQPSYPEQGY DRPYEDSSQH YYEGGNSQYG QQQDAYQGPP PQQGYPPQQQ QYPGQQGYPG QQQGYGPSQG GPGPQYPNYP QGQGQQYGGY RPTQPGPPQP PQQRPYGYDQ GQYGNYYQQ
Source	E.coli
Target Names	SS18
Protein Names	Recommended name: Protein SSXT Alternative name(s): Protein SYT Synovial sarcoma translocated to X chromosome protein
Expression Region	2-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.