



Recombinant Human Dixin (DIXDC1)

Product Code	CSB-EP621878HU-B
Abbreviation	DIXDC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q155Q3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLACLTRGNL LDVLQEGFNE QQLQAYVAWV NAQLKKRPAV KPVQDLRQDL RDGVILAYLI EIVAGEKLSG VQLSPGNQQE MKNNVEKVLQ FVASKKIRMH QTSAKDIVDG NLKSIMRLVL ALAAHFKPGS SRTVNQGRDS RAPLQSHRPH CATAVAQGAA AALADVCHDM SRSGRDVFRY RQRNSSMDEE IENPYWSVRA LVQQYEGQQR SPSESSCSSL TSPSPIHSAK SESIITQSEE KADFVIIPAE GIENRTEGTD SPLSRDWRPG SPGTYLETSW EEQLLEQQEY LEKEMEEAKK MISGLQALLL NGSLPEDEQE RPLALCEPGV NPPEQLIIQ SRLDQSMEEN QDLKKELLKC KQEARNLQGI KDALQQRLTQ QDTSVLQLKQ ELLRANMDKD ELHNQNVDLQ RKLDERNRLG GEYKKELGQK DRLLQQHQAK LEEALRKLSD VSYHQVDLER ELEHKDVLLA HCMKREADEA TNYNSHNSQS NGFLLPTAGK GATSVSNRGT SDLQLVRDAL RSLRNSFSGH DPQHHTIDSL EQGISSLMER LHVMEQKKQ ERKVRVKSPR TQVGSEYRES WPPNSKLPHS QSSPTVSSTC TKVLYFTDRS LTPFMVNIPK RLEEVTLKDF KAAIDREGNH RYHFKALDPE FGTVKEEIFH DDDAIPGWEG KIVAWVEEDH GEN
Source	E.coli
Target Names	DIXDC1
Protein Names	Recommended name: Dixin Alternative name(s): Coiled-coil protein DIX1 Short name= Coiled-coil-DIX1 DIX domain-containing protein 1
Expression Region	1-683
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.