



Recombinant Human Adipose most abundant gene transcript 2 protein (APM2)

Product Code	CSB-BP621887HU
Abbreviation	APM2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15847
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASKGLQDLK QQVEGTAQEA VSAAGAAAQQ VVDQATEAGQ KAMDQLAKTT QETIDKTANQ ASDTFSGIGK KFGLLK
Source	Baculovirus
Target Names	ADIRF
Protein Names	Recommended name: Adipose most abundant gene transcript 2 protein Short name= apM-2
Expression Region	1-76
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	APM2 gene is exclusively expressed in adipose tissue. Its function is currently unknown.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.