



Recombinant Human Zinc finger and BTB domain-containing protein 6 (ZBTB6)

Product Code	CSB-MP621889HU
Abbreviation	ZBTB6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15916
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAESDVLHF QFEQQGDVVL QKMNLLRQQN LFCDVSIYIN DTEFQGHKVI LAACSTFMRD QFLLTQSKHV RITILQSAEV GRKLLLSCYT GALEVKRKEL LKYLTAASYL QMVHIVEKCT EALSKYLEID LSMKNNNQHT DLCQSSDPDV KNEDENSDKD CEIIEISEDS PVNIDFHVKE EESNALQSTV ESLTSERKEM KSPELSTVDI GFKDNEICIL HVESISTAGV ENGQFSQPCT SSKASMYFSE TQHSLINSTV ESRVAEVPNG QDQGLFCENT EGSYGTVSEI QNLEEGYSLR HQCPRCPRGF LHVENYLRHL KMHKLFLCLQ CGKTFTQKKN LNRHIRGHMG IRPFQCTVCL KTFTAKSTLQ DHLNIHSGDR PYKCHCCDMD FKHKSALKKH LTSVHGRSSG EKLSRPDLKR QSLI
Source	Mammalian cell
Target Names	ZBTB6
Protein Names	Recommended name: Zinc finger and BTB domain-containing protein 6 Alternative name(s): Zinc finger protein 482 Zinc finger protein with interaction domain
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.