



Recombinant Human Cleavage and polyadenylation specificity factor subunit 6 (CPSF6)

Product Code	CSB-YP621966HU
Abbreviation	CPSF6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16630
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MADGVDHIDI YADVGEEFNQ EAEYGGHDQI DLYDDVISPS ANNGDAPEDR DYMDTLPPTV GDDVGKGAAP NVVYTYTGKR IALYIGNLTW WTTDEDLTEA VHSLGVNDIL EIKFFENRAN GQSKGFALVG VGSEASSKKL MDLLPKRELH GQNPVVTPCN KQFLSQFEMQ SRKTTQSGQM SGEKAGPPG GSSRAAFPQG GRGRGRFPGA VPGGDRFPGP AGPGGPPPPF PAGQTPRPP LGPPGPPGPP GPPPPGQVLP PPLAGPPNRG DRPPPPVLF GQPFQPPPLG PLPPGPPPPV PGYGGPPGPP PPQQGPPPPP GPFPPRPPG LGPPLTLAPP PHLPGPPGA PPPAPHVNPA FFPPPTNSGM PTSDSRGPPP TDPYGRPPPY DRGDYGPGR EMDTARTPLS EAEFEEIMNR NRAISSAIS RAVSDASAGD YGSAIETLVT AISLIKQSKV SADDRCKVLI SSLQDCLHGI ESKSYGSGSR RERSRERDHS RSREKSRRHK SRSRDRHDDY YRERSRERER HRDRDRDRDR ERDREREYRH R
Source	Yeast
Target Names	CPSF6
Protein Names	Recommended name: Cleavage and polyadenylation specificity factor subunit 6 Alternative name(s): Cleavage and polyadenylation specificity factor 68 kDa subunit Short name= CFIm68 Short name= CPSF 68 kDa subunit Pre-mRNA cleavage
Expression Region	1-551
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is one subunit of a cleavage factor required for 3 RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of



the earliest steps in the assembly of the 3 end processing complex and facilitates the recruitment of other processing factors. The cleavage factor complex is composed of four polypeptides. This gene encodes the 68kD subunit. It has a domain organization reminiscent of spliceosomal proteins.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.