



Recombinant Human Nuclear factor interleukin-3-regulated protein (NFIL3)

Product Code	CSB-EP621967HU
Abbreviation	NFIL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16649
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MQLRKMQTVK KEQASLDASS NVDKMMVLNS ALTEVSEDST TGEELLLSEG SVGKNKSSAC RRKREFIPDE KKDAMYWEKR RKNNEAAKRS REKRRRLNDLV LENKLIALGE ENATLKAELL SLKLFGLIS STAYAQEIQK LSNSTAVYFQ DYQTSKSNVS SFVDEHEPSM VSSCISVIK HSPQSSLSDV SEVSSVEHTQ ESSVQGSCRS PENKFQIIKQ EPMELESYTR EPRDDRGSYT ASIQNYMGN SFGYSHSPP LLQVNRSSSN SPRTSETDDG VVGKSSDGED EQQVPKGPPIH SPVELKHVHA TVVKVPEVNS SALPHKLRIK AKAMQIKVEA FDNEFEATQK LSSPIDMTSK RHFELEKHS A PSMVHSSLTP FSVQVTNIQD WSLKSEHWHQ KELSGKTQNS FKTGVVEMKD SGYKVSDPEN LYLKQGIANL SAEVVSLKRL IATQPISASD SG
Source	E.coli
Target Names	NFIL3
Protein Names	Recommended name: Nuclear factor interleukin-3-regulated protein Alternative name(s): E4 promoter-binding protein 4 Interleukin-3 promoter transcriptional activator Interleukin-3-binding protein 1 Transcriptional activator NF-IL3A
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.