



Recombinant Uncharacterized protein B0035.15 (B0035.15)

Product Code	CSB-YP622014CXY
Abbreviation	B0035.15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17439
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MKRLWKKVGS GVKNANYLKK SVEEVAKIEK GETKVVPPKY PTEKSREVS EIKKELAMKN ESLVENMNKM YIKSTDPVER WTSTKDLPTR ESEFLHRNDP IWEYGFYEPP VERIPKDKLM FKEALEYLRQ RQELLSDSSS PAQKKQAAEF VADHVVTSRV NAETLDDIYE YFRPFERKDK QKVVNRHALA ALQDHSVQNS DERRILDEAK DVGKKIRHIS NNTKFLDEYQ RLEDEEKDKV REAIAQLRQE EHERLNKRLG QLGEIEKTAQ EAMKKAVEEK KEK
Source	Yeast
Target Names	B0035.15
Protein Names	Recommended name: Uncharacterized protein B0035.15
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.