



Recombinant Human Protein FAM98C (FAM98C)

Product Code	CSB-EP622041HU-B
Abbreviation	FAM98C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17RN3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEAVKAEawe GAAVAQDLLA LGYGGVPGAA SRGASCPDFR GLCVRLAAEL ATLGALEQQR EAGAEVLSAG DPGGAEEDFL RQLGSLREL HCPDRALCGG DGAAALREPG AGLRLLRFLC SELQATRLLC LRSLLDPSPR PPLGEGVVEG AGMVQELDLT LQALGLPRPA PGTPASQLLQ ELHAKISELQ PSLPPGSLQP LLSCSLDAPR WEALESLSQS LRDQYRCRRC LLLKRLDLTT SAFHWSdraE AQGEAMRAVL IPIREVLTPe SDISIAHVLA ARADLSCLVP ATSVAVRRGT CCAINKVLMG NVPDRGGRPN ELEPPMPTWR SRREDGGPQC WGRKKKKKK
Source	E.coli
Target Names	FAM98C
Protein Names	Recommended name: Protein FAM98C
Expression Region	1-349
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.