



Recombinant Human Zinc finger SWIM domain-containing protein 7 (ZSWIM7)

Product Code	CSB-EP622178HU
Abbreviation	ZSWIM7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19AV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAVVLPAVVE ELLSEMAAAV QESARIPDEY LLSLKFLFGS SATQALDLVD RQSITLISSP SGRRVYQVLG SSSKTYTCLA SCHYCSCPAF AFSVLRKSDS ILCKHLLAVY LSQVMRTCQQ LSVSDKQLTD ILLMEKKQEA
Source	E.coli
Target Names	ZSWIM7
Protein Names	Recommended name: Zinc finger SWIM domain-containing protein 7 Alternative name(s): SWIM domain-containing and Srs2-interacting protein 1 homolog
Expression Region	1-140
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.