



Recombinant Schizosaccharomyces pombe Serine/threonine-protein kinase crk1 (crk1)

Product Code	CSB-BP622439SXV
Abbreviation	crk1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12126
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MDIEKSDKWT YVKERKVGEG TYAVVFLGRQ KETNRRVAIK KIKVGQFKDG IDISALREIK FLRESRHDNV IELVDVFSTK SNLNIILEFL DSDLEMLIKD KFIVFQPAHI KSWMVMLLRG LHHIHSR FIL HRDLKPNLL ISSDGVLKLA DFGLSRDFGT PSHMSHQVIT RWYRPPPELFM GCRSYGTGVD MWSVGCIFAE LMLRTPYLPG ESDLQQLNVI FRALGTPEPE VIKSMQQLPN YVEMKHIPPP NGGMEALFSA AGHEEIDLLK MMLDYNPYRR PTAQQALEHH YFSALPKPTH PSSLPRKGG EGIKHVSSDL QRQNNFPMRA NIKFV
Source	Baculovirus
Target Names	crk1
Protein Names	Recommended name: Serine/threonine-protein kinase crk1 EC= 2.7.11.23 Alternative name(s): Mitotic catastrophe suppressor 6
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.