



# Recombinant *Saccharomyces cerevisiae* Transcriptional modulator WTM1 (WTM1)

<b>Product Code</b>	CSB-MP622467SVG
<b>Abbreviation</b>	WTM1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12363
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPKKVWKSST PSTYEHISL RPKFVSRVDN VLHQRKSLTF SNVVVPDKKN NLTSSVIYS QGSDIYEIDF AVPLQEAASE PVKDYGDAFE GIENTSLSPK FVYQGETVSK MAYLDKTGET TLLSMSKNGS LAWFKEGIKV PIHIVQELMG PATSYASIHSLTRPGDLPEK DFSLAISDFG ISNDTETIVK SQSNGDEEDS ILKIIDNAGK PGEILRTVHV PGTTVTHTVR FFDNHIFASC SDDNILRFWD TRTSDKPIWV LGEPKNGKLT SFDCSQVSNN LFVTGFSTGI IKLWDARAAE AATDLTYRQ NGEDPIQNEI ANFYHAGGDS VVDVQFSATS SSEFFTVGGT GNIYHWNTDY SLSKYNPDDT IAPPQDATEE SQTKSLRFLH KGGSRSPKQ IGRRNTAAWH PVIENTLVGTV DDDSLVSIYK PYTEESE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	WTM1
<b>Protein Names</b>	Recommended name: Transcriptional modulator WTM1
<b>Expression Region</b>	1-437
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.