



Recombinant *Saccharomyces cerevisiae* Serine/threonine-protein phosphatase 2A activator 2 (RRD2)

| | |
|--------------------------|---|
| Product Code | CSB-MP622483SVG |
| Abbreviation | RRD2 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q12461 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLPEKRLLLTP DDMKLWEESP TRAHFTKFII DLAESVKGHE NSQYKEPISE SINSMMNLLS QIKDITQKHP VIKDADSSRF GKVEFRDFYD EVSRNSRKIL RSEFPLTDE QLEQLSIYLD ESWGKRRID YGSGHELNFM CLLYGLYSYG IFNLSNDSTN LVLKVFIEYL KIMRILETKY WLEPAGSHGV WGLDDYHFLP FLFGAFQLTT HKHLKPISIH NNELVEMFAH RYLYFGCIAF INKVKSSASL RWHSPMLDDI SGVKTWSKVA EGMIKMYKAE VLSKLPIMQH FYFSEFLPCP DGVSPPRGHI HDGTDKDDEC NFEGHVHTTW GDCCGIKLPS AIAATEMNKK HHKPIPFD |
| Source | Mammalian cell |
| Target Names | RRD2 |
| Protein Names | Recommended name: Serine/threonine-protein phosphatase 2A activator 2 EC=5.2.1.8 Alternative name(s): Peptidyl-prolyl cis-trans isomerase PTPA-2 Short name= PPIase PTPA-2 Short name= Rotamase PTPA-2 Phosphotyrosyl phospho |
| Expression Region | 1-358 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.