



Recombinant *Saccharomyces cerevisiae* Protein ZPS1 (ZPS1)

Product Code	CSB-BP622490SVG
Abbreviation	ZPS1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12512
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	APVTYDTNST AELQSPSSQE ILGWSHATFP TIYQTCNETN ARMLNAAFKD TAEITAYGKD RLLNYGVDDV YYKRWFGNGS IFTVMGVFEQ LMEASKGAML MRCDDIDGLC AANPNYYAGH HRQSAPAETV ICDYFYTSKK PLSTICFEGT IVDVGPKHYA GIDMLHRYLH VPTMSMDGYV GEYAETLEEV VDYTQNNATY AVRNTDNYLY YLADVYSASV IPGGCLGNL
Source	Baculovirus
Target Names	ZPS1
Protein Names	Recommended name: Protein ZPS1
Expression Region	21-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.