



# Recombinant Human Semaphorin-3B (SEMA3B)

<b>Product Code</b>	CSB-MP622646HU
<b>Abbreviation</b>	SEMA3B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13214
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>APSPPR LRLSFQELQA WHGLQTFSL E RTCCYQALLV DEERGRLFVG  AENHVASLNL DNISKRAKKL AWPAPVEWRE ECNWAGKDIG TECMNFVKLL  HAYNRTHLLA CGTGAFHPTC AFVEVGHRAE EPVLRLLDPGR IEDGKKGKSPY  DPRHRAASVL VGEELYSGVA ADLMGRDFTI FRSLGQRPSL RTEPHDSRWL  NEPKFVKVFW IPESENPDDD KIYFFFRETA VEAAPALGRL SVSRVGGQICR  NDVGGQRSLV NKWTTFLKAR LVCSVPGVEG DTHFDQLQDV FLLSSRDHRT  PLLYAVFSTS SSIFQGSAVC VYSMNDVRRR FLGPF AHKEG PMHQWVSYQG  RVPYPRPGMC PSKTFGTFFS TKDFPDDVIQ FARNHPLMYN SVLPTGGRPL  FLQVGANYTF TQIAADRVA A ADGHYDVLFI GTDVGTVLKV ISVPGKSRPS  AEGLLLEELH VFEDSAAVTS MQISSKRHQL YVASRSAVAQ IALHRCAAHG  RVCTECCLAR DPYCAWDGVA CTRFQPSAKR RFRRQDVRNG  DPSTLCSGDS SRPALLEHKV FGVEGSSAFL ECEPRSLQAR VEWTFQRAGV  TAHTQVLAEE RTERTARGLL LRRLRRRDSG VYLCAAVEQG FTQPLRRLSL  HVLSATQAER LARAEAAAPA APPGPKLWYR DFLQLVEPPG GGSANSLRMC  RPQPALQSLP LESRRKGRNR RTHAPEPRAE RGPRSATHW</p>
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SEMA3B
<b>Protein Names</b>	Recommended name: Semaphorin-3B Alternative name(s): Sema A(V) Semaphorin-V Short name= Sema V
<b>Expression Region</b>	25-749
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	The semaphorin/collapsin family of molecules plays a critical role in the guidance of growth cones during neuronal development. The secreted protein encoded by this gene family member is important in axonal guidance and has been shown to act as a tumor suppressor by inducing apoptosis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.