



Recombinant Human Ecto-ADP-ribosyltransferase 3 (ART3)

Product Code	CSB-EP622665HU-B
Abbreviation	ART3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13508
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	EVLD MADNAFDDEY LKCTDRMEIK YVPQLLKEEK ASHQQLDTVW ENAKAKWAAR KTQIFLPMNF KDNHGIALMA YISEAQEQTP FYHLFSEAVK MAGQSREDYI YGFQFKAFHF YLTRALQLLR KPCEASSKTV VYRTSQGTSF TFGGLNQARF GHFTLAYS AK PQAANDQLTV LSIYTCLGVD IENFLDKESE RITLIPLNEV FQVSQEGAGN NLILQSINKT CSHYECAFLG GLKTENCIEN LEYFQPIYVY NPGEKNQKLE DHSEKNWKLE DHGEKNQKLE DHGVKILEPT QIPGMKIPEP FPLPEDKSQG NINNPTPGPV PVPGPKSHPS AS
Source	E.coli
Target Names	ART3
Protein Names	Recommended name: Ecto-ADP-ribosyltransferase 3 EC= 2.4.2.31 Alternative name(s): Mono(ADP-ribosyl)transferase 3 NAD(P)(+)-arginine ADP-ribosyltransferase 3
Expression Region	27-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.