



Recombinant Human BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 (TNFAIP1)

Product Code	CSB-MP622686HU
Abbreviation	TNFAIP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13829
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSGDTCLCPA SGAKPKLSGF KGGGLGNKYV QLNVGGSLLYY TTVRALTRHD TMLKAMFSGR MEVLTDKEGW ILIDRCGKHF GTILNYLRDD TITLPQNRQE IKELMAEAKY YLIQGLVNM C QSALQDKKDS YQPVCNIPII TSLKEEERLI ESSTKPVVKL LYNRSNNKYS YTSNSDDHLL KNIELFDKLS LRFNGRVLFI KDVIGDEICC WSFYGQGRKL AEOVCTSIVY ATEKKQTKVE FPEARIYEET LNVLLYETPR VPDNSLLEAT SRSRSQASPS EDEETFELRD RVRRIHVKRY STYDDRQLGH QSTHRD
Source	Mammalian cell
Target Names	TNFAIP1
Protein Names	Recommended name: BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 Short name= hBACURD2 Alternative name(s): BTB/POZ domain-containing protein TNFAIP1 Protein B12 Tumor necrosis factor, alpha-induced pr
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene was identified as a gene whose expression can be induced by the tumor necrosis factor alpha (TNF) in umbilical vein endothelial cells. Studies of a similar gene in mouse suggest that the expression of this gene is developmentally regulated in a tissue-specific manner.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.