



# Recombinant Human BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 (TNFAIP1)

<b>Product Code</b>	CSB-MP622686HU
<b>Abbreviation</b>	TNFAIP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q13829
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSGDTCLCPA SGAKPKLSGF KGGGLGNKYV QLNVGGSLLYY TTVRALTRHD TMLKAMFSGR MEVLTDKEGW ILIDRCGKHF GTILNYLRDD TITLPQNRQE IKELMAEAKY YLIQGLVNM C QSALQDKKDS YQPVCNIPII TSLKEEERLI ESSTKPVVKL LYNRSNNKYS YTSNSDDHLL KNIELFDKLS LRFNGRVLFI KDVGDEICC WSFYGQGRKL AEOVCTSIVY ATEKKQTKVE FPEARIYEET LNVLLYETPR VPDNSLLEAT SRSRSQASPS EDEETFELRD RVRRIHVKRY STYDDRQLGH QSTHRD
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TNFAIP1
<b>Protein Names</b>	Recommended name: BTB/POZ domain-containing adapter for CUL3-mediated RhoA degradation protein 2 Short name= hBACURD2 Alternative name(s): BTB/POZ domain-containing protein TNFAIP1 Protein B12 Tumor necrosis factor, alpha-induced pr
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene was identified as a gene whose expression can be induced by the tumor necrosis factor alpha (TNF) in umbilical vein endothelial cells. Studies of a similar gene in mouse suggest that the expression of this gene is developmentally regulated in a tissue-specific manner.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.