



Recombinant Human Krueppel-like factor 9 (KLF9)

| | |
|--------------------------|---|
| Product Code | CSB-BP622687HU |
| Abbreviation | KLF9 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q13886 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSAAAYMDFV AAQCLVSISN RAAVPEHGVA PDAERLRLPE REVTKEHGDP GDTWKDYCTL VTIAKSLLDL NKYRPIQTPS VCSDSLESPD EDMGSDSDVT TESGSSPSHS PEERQDPGSA PSPLSLLHPG VAAKGKHASE KRHKCPYSGC GKVYGKSSHL KAHYRVHTGE RPFPCWPDC LKKFSRSDEL TRHYRHTTGE KQFRCPLCEK RFMRSDHLTK HARRHTEFHP SMIKRSKKAL ANAL |
| Source | Baculovirus |
| Target Names | KLF9 |
| Protein Names | Recommended name: Krueppel-like factor 9 Alternative name(s): Basic transcription element-binding protein 1 Short name= BTE-binding protein 1 GC-box-binding protein 1 Transcription factor BTEB1 |
| Expression Region | 1-244 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Target Details | This protein is a transcription factor that binds to GC box elements located in the promoter. Binding of the encoded protein to a single GC box inhibits mRNA expression while binding to tandemly repeated GC box elements activates transcription. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |