



Recombinant Human Conserved oligomeric Golgi complex subunit 2 (COG2)

Product Code	CSB-YP622787HU
Abbreviation	COG2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14746
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEKSRMNLPK GPDTLCFDKD EFMKEDFDVD HFVSDCRKRV QLEELRDDLE LYYKLLKTAM VELINKDYAD FVNLSTNLVG MDKALNQLSV PLGQLREEVL SLRSSVSEGI RAVDERMSKQ EDIRKKKMCV LRLIQVIRSV EKIEKILNSQ SSKETSALEA SSPLLTGQIL ERIATEFNQL QFHAVQSKGM PLLDKVRPRI AGITAMLQQS LEGLLLEGLQ TSDVDIIRHC LRTYATIDKT RDAEALVGQV LVKPYIDEVI IEQFVESHVN GLQVMYNKLL EFVPHHCRLL REVTGGAISS EKGNTVPGYD FLVNSVWPQI VQGLEEKLPS LFNPGNPDAF HEKYTISMDF VRRLERQCGS QASVKRLRAH PAYHSFNKKW NLPVYFQIRF REIAGSLEAA LTDVLEDAPA ESPYCLLASH RTWSSLRRCW SDEMFLPLL V HRLWRLTLQI LARYSVFVNE LSLRPISNES PKEIKKPLVT GSKEPSITQG NTEDQGSGPS ETKPVVSISR TQLVYVVADL DKLQEQLPEL LEIIPKLEM IGFKNFSSIS AALEDSQSSF SACVPSLSSK IIQDLSDSCF GFLKSALEVP RLYRRTNKEV PTTASSYVDS ALKPLFQLQS GHKDKLKQAI IQQWLEGTLS ESTHKYYETV SDVLNSVKKM EESLKRKQA RKTTPANPVG PSGGMSDDDK IRLQLALDVE YLGEQIQKLG LQASDIKSFS ALAELVAAAK DQATAEQP
Source	Yeast
Target Names	COG2
Protein Names	Recommended name: Conserved oligomeric Golgi complex subunit 2 Short name= COG complex subunit 2 Alternative name(s): Component of oligomeric Golgi complex 2 Low density lipoprotein receptor defect C-complementing protein
Expression Region	1-738
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.