



Recombinant Human Coiled-coil domain-containing protein 6 (CCDC6)

Product Code	CSB-YP623068HU
Abbreviation	CCDC6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16204
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> ADSASESDT DGAGGNSSSS AAMQSSCSST SGGGGGGGGG GGGGKSGGIV ISPFRLEELT NRLASLQQEN KVLKIELETY KLKCKALQEE NRDLRKASVT IQARAEQEEE FISNTLFKKI QALQKEKETL AVNYEKEEEF LTNELSRKLM QLQHEKAELE QHLEQEDEFQ VNKLMKKIKK LENDTISKQL TLEQLRREKI DLENTLEQEQ EALVNRLWKR MDKLEAEKRI LQEKLDQPVS APPSPRDISM EIDSPENMMR HIRFLKNEVE RLKKQLRAAQ LQHSEKMAQY LEEERHMREE NLRLQRKLQR EMERREALCR QLSESESSLE MDDERYFNEM SAQGLRPRTV SSSIPYTPSP SSSRPISPGL SYASHTVGFT PPTSLTRAGM SYYNSPGLHV QHMGTSHGIT RPSPRRSNSP DKFKRTPPPP SPNTQTPVQP PPPPPPPMQ PTVPSAATSQ PTPSQHSAHP SSQP </p>
Source	Yeast
Target Names	CCDC6
Protein Names	Recommended name: Coiled-coil domain-containing protein 6 Alternative name(s): Papillary thyroid carcinoma-encoded protein Protein H4
Expression Region	2-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.