



# Recombinant Human snRNA-activating protein complex subunit 1 (SNAPC1)

<b>Product Code</b>	CSB-EP623083HU
<b>Abbreviation</b>	SNAPC1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16533
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGTPPGLQTD CEALLSRFQE TDSVRFEDFT ELWRNMKFGT IFCGRMRNLE KNMFTKEALA LAWRYFLPPY TFQIRVGALY LLYGLYNTQL CQPKQKIRVA LKDWDEVLFK QDDLVAQHF DAAYIFRKLK LDRAFHFTAM PKLLSYRMKK KIHRAEVTEE FKDPDRVMK LITSDVLEEM LNVHDHYQNM KHVISVDKSK PDKALSLIKD DFFDNIKNIV LEHQQWHKDR KNPSLKSSTN DGEEKMEGNS QETERCERAE SLAKIKSKAF SVVIQASKSR RHRQVKLDSS DSDSASGGQQ VKATRKKEKK ERLKPAGRKM SLRNKGNVQN IHKEDKPLSL SMPVITEEEE NESLSGTEFT ASKKRRKH
<b>Source</b>	E.coli
<b>Target Names</b>	SNAPC1
<b>Protein Names</b>	Recommended name: snRNA-activating protein complex subunit 1 Short name= SNAPc subunit 1 Alternative name(s): Proximal sequence element-binding transcription factor subunit gamma Short name= PSE-binding factor subunit gamma Short
<b>Expression Region</b>	1-368
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.