



Recombinant Human snRNA-activating protein complex subunit 1 (SNAPC1)

Product Code	CSB-MP623083HU
Abbreviation	SNAPC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16533
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGTPPGLQTD CEALLSRFQE TDSVRFEDFT ELWRNMKFGT IFCGRMRNLE KNMFTKEALA LAWRYFLPPY TFQIRVGALY LLYGLYNTQL CQPKQKIRVA LKDWDDEVLFK QQDLVNAQHF DAAYIFRKLK LDRAFHFTAM PKLLSYRMKK KIHRAEVTEE FKDPSTRVMK LITSDVLEEM LNVHDHYQNM KHVISVDKSK PDKALSLIKD DFFDNIKNIV LEHQQWHKDR KNPSLKSSTN DGEEKMEGNS QETERCERAE SLAKIKSKAF SVVIQASKSR RHRQVKLDSS DSDSASGGQQ VKATRKKEKK ERLKPAGRKM SLRNKGNVQN IHKEDKPLSL SMPVITEEEE NESLSGTEFT ASKKRRKH
Source	Mammalian cell
Target Names	SNAPC1
Protein Names	Recommended name: snRNA-activating protein complex subunit 1 Short name= SNAPc subunit 1 Alternative name(s): Proximal sequence element-binding transcription factor subunit gamma Short name= PSE-binding factor subunit gamma Short
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.