



Recombinant Human Serotonin N-acetyltransferase (AANAT)

Product Code	CSB-MP623090HU
Abbreviation	AANAT
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16613
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSTQSTHPLK PEAPRLPPGI PESPSCQRRH TLPASEFRCL TPEDAVSAFE IEREAFISVL GVCPLYLDEI RHFLTLCPEL SLGWFEEGCL VAFIIGSLWD KERLMQESLT LHRSGGHIAH LHVLAHVRAF RQQGRGPILL WRYLHHLGSQ PAVRRALMC EDALVPFYER FSFHAVGPCA ITVGSLTFME LHCSLRGHPF LRRNSGC
Source	Mammalian cell
Target Names	AANAT
Protein Names	Recommended name: Serotonin N-acetyltransferase Short name= Serotonin acetylase EC= 2.3.1.87 Alternative name(s): Aralkylamine N-acetyltransferase Short name= AA-NAT
Expression Region	1-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein belongs to the acetyltransferase superfamily. It is the penultimate enzyme in melatonin synthesis and controls the night/day rhythm in melatonin production in the vertebrate pineal gland. Melatonin is essential for the function of the circadian clock that influences activity and sleep. This enzyme is regulated by cAMP-dependent phosphorylation that promotes its interaction with 14-3-3 proteins and thus protects the enzyme against proteasomal degradation. This gene may contribute to numerous genetic diseases such as delayed sleep phase syndrome. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.