



Recombinant WD repeat and FYVE domain-containing protein 2 (wdfy-2)

Product Code	CSB-EP623235CXY-B
Abbreviation	wdfy-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q18964
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAAIINQRVE QGESSMGGAK PALLHRITGH VARINDVILL SKDEGVWTAS DDRSVRLYLK RDNDQFWPSI HHFMPVAPTS LYYSEETYKL LVGLINGNVY EFSVADDFNS MTESRKWTCH AGPISGLGFA LSSELIFSCS RDKSIVWHCS ENSNKMGSYL LENSCTAMVI DLPFVFGVDH GGHVTVLRIV DNQPNLVSKL SAHTNAITSL TWDGNKKVLY SGSSDHLIIM WDIGGGRGEA YELNGHNGKV TTLCAAPAAK RLFSADDEHGK LMCWDMNCKR VETPEWKTSD CCQKCNQPPF WNLQAMWQRK VVGLRQHHCR TCGSAVCGSC CDNWTTYPPM GYETKIRICN DCNARMKDNP QNFNLTP LAI PHEIHTGITA MHLQETLGLL VTSGQNRVIM IWDVRSVCSA PSGSQ
Source	E.coli
Target Names	wdfy-2
Protein Names	Recommended name: WD repeat and FYVE domain-containing protein 2
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.