



Recombinant *Saccharomyces cerevisiae* Histone deacetylase HOS1 (HOS1)

| | |
|--------------------------|---|
| Product Code | CSB-MP623590SVG |
| Abbreviation | HOS1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q12214 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSKLVISTSI FQSQVADLLP CNNHQKSQLT YSLINAYDLL QHFDEVLTFF YARKDDLLEF HSKSYIDYLI NGRFNKMMQA DVNNPMVESK WSELSELADN WNEKIDYNPS QDLQRFTTRE NLYNYLNS QALENNMDCI NNSEVPTNDK PTDTYILNSE TKQYNLEGDC PIFSYPMYC QVITGATLNL LDHLSPTERL IGINWDGGRH HAFKQRASGF CYINDVLLI QRLRKAKLNK ITYVDFDLHH GDGVEKAFQY SKQIQTISVH LYEPGFFPGT GSLSDSRKDK NVVNIPLKHG CDDNYLELIA SKIVNPLIER HEPEALIEC GGDGLLGDRF NEWQLTIRGL SRIINIMKS YPRAHIFLLG GGGYNDLLMS RFYTYLTWCV TKQFSNLRCG DNNSFQIDPF DVCDGDDSEQ FIREHDLVEM YNEENYQYWI YEMEGSSRMK MLRNDNKDRD MVELMKFYEL |
| Source | Mammalian cell |
| Target Names | HOS1 |
| Protein Names | Recommended name: Histone deacetylase HOS1 EC= 3.5.1.98 |
| Expression Region | 1-470 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.