



Recombinant Endo-1,4-beta-xylanase 4 (XYN4)

Product Code	CSB-BP623632DOZ
Abbreviation	XYN4
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q12550
Product Type	Recombinant Protein
Immunogen Species	Aspergillus niger
Purity	>85% (SDS-PAGE)
Sequence	PPVP EPVLVRSRAG INYVQNYNGN LGDFTYDESA GTFSMYWEDG VSSDFVVGLG WTTGSSNAIT YSAEYSASGS SSYLAVYGWV NLSQAEYYIV EDYGDYNPCS SATSLGTEYS DGSTYQVCTD TRTNEPSITG TSTFTQYFSV RESTRTSGTV TVAIHFNFVA QHGFNSDFN YQVMAVEAWS GACSASVTIS S
Source	Baculovirus
Target Names	XYN4
Protein Names	Recommended name: Endo-1,4-beta-xylanase 4 Short name= Xylanase 4 EC= 3.2.1.8 Alternative name(s): 1,4-beta-D-xylan xylanohydrolase 4
Expression Region	17-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>