



Recombinant Human Secernin-1 (SCRN1)

Product Code	CSB-EP623649HU-B
Abbreviation	SCRN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12765
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAPPSYCF VAFPPRAKDG LVVFGKNSAR PRDEVQEVVY FSAADHEPES KVECTYISID QVPRTYAIMI SRPAWLWGAE MGANEHGVC I ANEAINREP AAEIEALLGM DLVRLGLERG ETAKEALDVI VSLLEEHGQG GNYFEDANSC HSFQSAYLIV DRDEAWVLET IGKYWAAEKV TEGVRCICSQ LSLTTKMDAE HPELSYAQS QGWWTGE GEF NFSEVFSPVE DHLDCGAGKD SLEKQEE SIT VQTMNTLRD KASGVCIDSE FFLTTASGVS VLPQNRSSPC IHYFTGTPDP SRSIFKPFIF VDDVKLVPKT QSPCFGDDDD AKKEPRFQEK PDRRHELYKA HEWARAIIES DQEQRKLR S TMLELEKQGL EAMEEILTSS EPLDPAEVGD LFYDCVDTEI KFFK
Source	E.coli
Target Names	SCRN1
Protein Names	Recommended name: Secernin-1
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene likely encodes a member of the secernin family of proteins. A similar protein in rat functions in regulation of exocytosis in mast cells. Alternatively spliced transcript variants have been described.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.