




Recombinant Human Squamous cell carcinoma antigen recognized by T-cells 3 (SART3), partial

Product Code	CSB-EP624010HU
Relevance	U6 snRNP-binding protein that functions as a recycling factor of the splicing machinery. Promotes the initial reassembly of U4 and U6 snRNPs following their ejection from the spliceosome during its maturation. Also binds U6atac snRNPs and may function as a recycling factor for U4atac/U6atac spliceosomal snRNP, an initial step in the assembly of U12-type spliceosomal complex. The U12-type spliceosomal complex plays a role in the splicing of introns with non-canonical splice sites. May also function as a substrate-targeting factor for deubiquitinases like USP4 and USP15. Recruits USP4 to ubiquitinated PRPF3 within the U4/U5/U6 tri-snRNP complex, promoting PRPF3 deubiquitination and thereby regulating the spliceosome U4/U5/U6 tri-snRNP spliceosomal complex disassembly. May also recruit the deubiquitinase USP15 to histone H2B and mediate histone deubiquitination, thereby regulating gene expression and/or DNA repair. May play a role in hematopoiesis probably through transcription regulation of specific genes including MYC.
Abbreviation	Recombinant Human SART3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15020
Alias	Tat-interacting protein of 110 kDa1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QRKRARAEEKKALKKKKKKIRGPEKRGADEDDEKEWGDDEEEQPSKRRRVENSI PAAGETQNVEVAAGPAGKCAAVDVEPPSKQKEKAASLKRDMPPKVLHDSSKDS ITVFVSNLPYSMQEPDTKLRLPLFEACGEVVQIRPIFSNRGDFRGYCYVEFKEEK SALQALEMDRKSVEGRPMFVSPCVDKSKNPDFKVFYRSTSLKHKLFISGLPF SCTKEELEEEICKAHGTVKDLRLVTNRAGKPKGLAYVEYENESQASQAVMKMD GMTIKENIIKVAISNPPQRKVPEKPETRKAPGGPMLLP
Research Area	Cancer
Source	E.coli
Target Names	SART3
Protein Names	Recommended name: Squamous cell carcinoma antigen recognized by T-cells 3 Short name= SART-3 Short name= hSART-3 Alternative name(s): Tat-interacting protein of 110 kDa Short name= Tip110




Expression Region	600-900aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	49.8kDa
Protein Length	Partial


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP624010HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) SART3.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP624010HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) SART3.

Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.