



# Recombinant Human Cdc42-interacting protein 4 (TRIP10)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP624029HU-B   |
| <b>Abbreviation</b>      | TRIP10   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q15642   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | MDWGTTELWDQ FEVLERHTQW GLDLLDRYVK FVKERTEVEQ AYAKQLRSLV<br>KKYLPKRPAK DDPEKFSQQ QSFVQILQEV NDFAGQRELV AENLSVRVCL<br>ELTKYSQEMK QERKMHFQEG RRAQQQLENG FKQLENSKRK<br>FERDCREAEK AAQTAERLDQ DINATKADVE KAKQQAHLRS HMAEESKNEY<br>AAQLQRFNRD QAHFYFSQMP QIFDKLQDMD ERRATRLGAG YGLLSEAELE<br>VVPPIAKCLE GMKVAANAVID PKNDSHVLIE LHKSGFARPG DVEFEDFSQP<br>MNRAPSDSSL GTPSDGRPEL RGPGRSRTKR WPFGKKNKPR<br>PPPLSPLGGP VPSALPNGPP SPRSGRDPLA ILSEISKSVK PRLASFRSLR<br>GSRGTVVTEG FSHLPPEQQR KRLQQQLEER SRELQKEVDQ REALKMKMDV<br>YEKTPQMGGP ASLEPQIAET LSNIERLKLE VQKYEAWLAE AESRVLSNRG<br>DSLRSRHARPP DPPASAPPDS SSNSASQDTK ESSEPPSEE SQDTPIYTEF<br>DEDFEEPTS PIGHCVAIYH FEGSSEGTIS MAEGEDLSLM EEDKGDGWTR<br>VRRKEGGEGY VPTSYL RVTL N |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | TRIP10   |
| <b>Protein Names</b>     | Recommended name: Cdc42-interacting protein 4 Alternative name(s): Protein Felic Salt tolerant protein Short name= hSTP Thyroid receptor-interacting protein 10 Short name= TR-interacting protein 10 Short name= TRIP-10  |
| <b>Expression Region</b> | 1-601  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.