



Recombinant Human Cdc42-interacting protein 4 (TRIP10)

Product Code	CSB-YP624029HU
Abbreviation	TRIP10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15642
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDWGTSLWDQ FEVLERHTQW GLDLLDRYVK FVKERTEVEQ AYAKQLRSLV KKYLPKRPAK DDPESKFSQQ QSFVQILQEV NDFAGQRELV AENLSVRVCL ELTKYSQEMK QERKMHFQEG RRAQQQLENG FKQLENSKRK FERDCREAEK AAQTAERLDQ DINATKADVE KAKQQAHLRS HMAEESKNEY AAQLQRFNRD QAHFYFSQMP QIFDKLQDMD ERRATRLGAG YGLLSEAELE VVPPIAKCLE GMKVAANAVD PKNDSHVLIE LHKSGFARPG DVEFEDFSQP MNRAPSDSSL GTPSDGRPEL RGPGRSRTKR WPFGKKNKPR PPPLSPLGGP VPSALPNGPP SPRSGRDPLA ILSEISKSVK PRLASFRLSLR GSRGTVVTED FSHLPPEQQR KRLQQQLEER SRELQKEVDQ REALKKMKDV YEKTPQMGGP ASLEPQIAET LSNIERLKLE VQKYEAWLAE AESRVLSNRG DSLRSRHARPP DPPASAPPDS SSNSASQDTK ESSEPPSEE SQDTPIYTEF DEDFEEPTS PIGHCVAIYH FEGSSEGTIS MAEGEDLSLM EEDKGDGWTR VRRKEGGEGY VPTSYLRVTL N
Source	Yeast
Target Names	TRIP10
Protein Names	Recommended name: Cdc42-interacting protein 4 Alternative name(s): Protein Felic Salt tolerant protein Short name= hSTP Thyroid receptor-interacting protein 10 Short name= TR-interacting protein 10 Short name= TRIP-10
Expression Region	1-601
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.