



Recombinant Human ALX homeobox protein 1 (ALX1)

Product Code	CSB-MP624031HU
Abbreviation	ALX1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15699
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEFLSEKFAL KSPPSKNSDF YMGAGGGLEH VMETLDNESF YSKASAGKCV QAFGPLPRAE HHVRLERTSP CQDSSVNYGI TKVEGQPLHT ELNRAMDNCN SLRMSPVKGM QEKGELDELG DKCDSNVSSS KKRRHRTTFT SLQLEELEKV FQKTHYPDVY VREQLALRTE LTEARVQVWF QNRRAKWRKR ERYGQIQQAK SHFAATYDIS VLPRTDSYPQ IQNNLWAGNA SGGSVVTSCM LPRDTSSCMT PYSHSPRTDS SYTGFSNHQN QFSHVPLNNF FTDSLTLTGAT NGHAFETKPE FERRSSSIIV LRMKAKEHTA NISWAM
Source	Mammalian cell
Target Names	ALX1
Protein Names	Recommended name: ALX homeobox protein 1 Alternative name(s): Cartilage homeoprotein 1 Short name= CART-1
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The specific function of this gene has yet to be determined in humans; however, in rodents, it is necessary for survival of the forebrain mesenchyme and may also be involved in development of the cervix. Mutations in the mouse gene lead to neural tube defects such as acrania and meroanencephaly.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.