



# Recombinant Human Hydroxyacylglutathione hydrolase, mitochondrial (HAGH)

<b>Product Code</b>	CSB-YP624118HU
<b>Abbreviation</b>	HAGH
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16775
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AALGAAC ARRGLGPALL GVFCHTDLRK NLTVDEGTMK VEVLPALTDN YMYLVIDDET KEAAIVDPVQ PQKVVDAARK HGVKLTTVLT THHHWDHAGG NEKLVKLESG LKVYGGDDRI GALTHKITHL STLQVGSLSNV KCLATPCHTS GHICYFVSKP GGSEPPAVFT GDTLFVAGCG KFYEGTADM CKALLEVLGR LPPDTRVYCG HEYTINNLKF ARHVEPGNAA IREKLAWAKE KYSIGEPTVP STLAEFTYN PFMRVREKTV QQHAGETDPV TTMRAVRREK DQFKMPRD
<b>Source</b>	Yeast
<b>Target Names</b>	HAGH
<b>Protein Names</b>	Recommended name: Hydroxyacylglutathione hydrolase, mitochondrial EC= 3.1.2.6 Alternative name(s): Glyoxalase II Short name= Glx II
<b>Expression Region</b>	14-308
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	The enzyme encoded by this gene is classified as a thiolesterase and is responsible for the hydrolysis of S-lactoyl-glutathione to reduced glutathione and D-lactate.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.