



Recombinant Human Tumor protein D53 (TPD52L1)

Product Code	CSB-MP624123HU
Abbreviation	TPD52L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16890
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEAAQAGLLETEPLQGTDEDAVASADFSSMLSEEEKEELKAELVQLEDEITTL RQVLSAKERHLVEIKQKLGMLMNELKQNFSSKSWHDMQTTTAYKKHETLSH AGQKATAAFSNVGTAKKFGDMSYSIRHSISMPAMRNSPTFKSFEERVETTV TSLKTKVGGTNPNGGSFEEVLSSTAHASAQSLAGGSRRTKEEELQC
Source	Mammalian cell
Target Names	TPD52L1
Protein Names	Recommended name: Tumor protein D53 Short name= hD53Alternative name(s): Tumor protein D52-like 1
Expression Region	1-144aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Isoform 2
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.