



Recombinant Uncharacterized protein B0564.3 (B0564.3)

Product Code	CSB-EP624164CXY
Abbreviation	B0564.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17528
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MTINYHKEIM TSHPWTFLL LFKWKGSIWK AVYMETIIFL ICYGIISVIY KTAMGESSQR TFESLVRYFD KRLSYIPLEF VLGFFVTTVV NRWTKLYQTI GFIDNVGLMA NCYIRGATEK ARIYRRNIMR YCELVQILVF RDMSMRTRRR FPTMETVVAA GFMNKHELEL YNSYDTKYNS KLGTKYWIPA NWALCMTYKA RKDGYIESDY FKAQMEGEIR TWRTNIEWVC NYDWVPLPLM YPQLVCLAVN LYFLVSIAR QLVIEKHKMV DEVDVYFPVM TFLQFIFYMG WLKVIDVMLN PFGEDDDDFE TNALIDRNIT MGLMIADNPM STELPRKDPF YDEVDVPLLY SEESSNIPNH HYHGSVSEVR LEQKGNAPVM MMPHSQSAAN LRRMMSFKSV DEDEKDINAF SMSHDDARMR NWREVSLDSS FLADLNENKE WKIPTNPQKF
Source	E.coli
Target Names	best-1
Protein Names	Recommended name: Uncharacterized protein B0564.3
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.