



Recombinant Putative aspartyl aminopeptidase (F01F1.9)

Product Code	CSB-BP624304CXY
Abbreviation	F01F1.9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19087
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MAAALKPSAP EIRKAAQEFI NYLNKAVTPF HATQEVKDRL LQAGFTELPE SGHWDIQPTS KYFVTKNRS ILAFVGGSY KPGSGFSIVV GHTDSPCLRV KPISHQKSDK FLQVGVSTYG GGIWRTWFDR DLSVAGLVIV KNGEKLQHKL IDVKKPVLFI PNLAIHLETD RTTFKPNTET ELRPILETF AAGINAPQKP ESTGFADPRN ITNNHHPQFL GLIAKEAGCQ PEDIVDL DLY LYDTNKAAIV GMEDEFISGA RLDNQVGTYT AISGLLESLT GESFKNDPQI RIAACFDNEE VGSDSANGAS SSFTEFVLR LRSAGGSTTAF EEAIGKSM LI SADQAHATHP NYS AKHEENH RPAFHGGVVV KVN VNQRYAT TSTTHAALKQ VAFE AQVPLQ VVVVRNDSPC GSTVGPILAT KLGLQTV DVG CPQLAMHSIR EFADTSSIIYQ ATTLYSTFYE RLSTVLSNMQ
Source	Baculovirus
Target Names	dnpp-1
Protein Names	Recommended name: Putative aspartyl aminopeptidase EC= 3.4.11.21
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.