



Recombinant Cercopithecus pygerythrus Forkhead box protein G1 (FOXG1)

Product Code	CSB-MP624340CGX
Abbreviation	FOXG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1A1A5
Product Type	Recombinant Protein
Immunogen Species	Chlorocebus pygerythrus (Vervet monkey) (Cercopithecus pygerythrus)
Purity	>85% (SDS-PAGE)
Sequence	MLDMGDRKEV KMIPKSSFSI NSLVPEAVQN DNHHASHGHH NSHHPQHSHHH HHHHHHHPPP PAPQPPPPPPQ QQPPPPPPPP APQPPQTRGA PAADDDKGPQ QLLLPPPPPP PPAAALDGAK ADGLGGKGEP GGGPGELAPV GPDEKEKGAG AGGEEKKGAG EGGKDGEGGK EGEKKNKYE KPPFSYNALI MMAIRQSPEK RLTLNGIYEF IMKNFPYYRE NKQGWQNSIR HNLSLNKCFV KVPRHYDDPG KGNYWMLDPS SDDVFIGGTT GKLRRRSTTS RAKLAFKRGAL RTSTGLTFM DRAGSLYWPM SPFLSLHHPR ASSTLSYNGT TSAYPSHPMP YSSVLTQNSL GNNHSFSTAN GLSVDRLVNG EIPYATHHLT AAALAASVPC GLSVPCSGTY SLNPCSVNLL AGQTSYFFPH VPHPSMSTSQS STSMSARAAS SSTSPQAPST LPCESLRPSL PSFTTGLSGG LSDYFTHQNQ GSSSNPLIH
Source	Mammalian cell
Target Names	FOXG1
Protein Names	Recommended name: Forkhead box protein G1 Short name= FoxG1
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.