



Recombinant Danio rerio Tumor necrosis factor, alpha-induced protein 8-like protein 2 A (tnfaip8l2a)

Product Code	CSB-BP624751DIL
Abbreviation	tnfaip8l2a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1ECV8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MEAFSSKDMA MQAQKKILSR MASKSVAQMF IDDNSSEVLD ELYRVSKEHS GNRAEAQKV VKNMIKIAVKV GVLFRHEKFS AEELSIAQEF RKKLHHGAMT AISFQEVEFT FDKSVMINLL SDCRDLLLKL VEKHLTPKSL DRIRHVFNHY SNPELLTHLY DPQGNLWSNL SKICSGLNRM IEEGKL
Source	Baculovirus
Target Names	tnfaip8l2a
Protein Names	Recommended name: Tumor necrosis factor, alpha-induced protein 8-like protein 2 A Short name= TIPE2 A Short name= TNF alpha-induced protein 8-like protein 2 A Short name= TNFAIP8-like protein 2 A
Expression Region	1-186
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.