



# Recombinant Arabidopsis thaliana BTB/POZ and MATH domain-containing protein 5 (BPM5)

<b>Product Code</b>	CSB-MP625195DOA
<b>Abbreviation</b>	BPM5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1EBV6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSESVIQGSN PDRVLSPTSS KSVTQTVNGS HQFVIQGYSL AKGMGIGKHI ASDNFSVGGY QWGIFFYPDG KNPEDNSSYV SVFIALASEG TEVRALFELA LVDQSGKGGKH KVHSHFERSL DGGPYTLKYR GSMWGYKRFF RRSILETSDY LKDDCLIINC TVGVVVSEIL CPQLHSVHVP DSELGSHFGV LLDSMEGSDI TFNIAGEKFL AHKLVLAARS PFFKSKFFSE FEANNTTEVTI NDLEPKVFKA LLQFMYKDSL PEDVEPATAH TFERLKLSEI YETLIVKVLA AADKYDLIRL RLLCESHICK GSVKSVAKI LALADRYNAK ELKGVCLKFT AENLAAVLET DAYQQMKDEC VTLQSELLKA VAGHEEGSNS TGGAKSQSVW AQLSDGGGDT TSRHVRQRTT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	BPM5
<b>Protein Names</b>	Recommended name: BTB/POZ and MATH domain-containing protein 5 Alternative name(s): Protein BTB-POZ AND MATH DOMAIN 5 Short name= AtBPM5
<b>Expression Region</b>	1-410
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.